

产品名称: RAP1GAP Rabbit Monoclonal Antibody
产品货号: AMRe02522



产品概述 (Summary)

产品名称 (Production Name)	RAP1GAP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAP1GAP
别名 (Alternative Names)	RAP1GAP; KIAA0474; RAP1GA1; Rap1 GTPase-activating protein 1; Rap1GAP; Rap1GAP1
基因 ID (Gene ID)	5909
蛋白 ID (SwissProt ID)	P47736.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 73 kDa; Observed MW: 95 kDa

研究背景 (Background)

产品名称: RAP1GAP Rabbit Monoclonal Antibody
产品货号: AMRe02522




GTPase activator for the nuclear Ras-related regulatory protein RAP-1A (KREV-1), converting it to the putatively inactive GDP-bound state.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of RAP1GAP in mouse brain lysates using RAP1GAP antibody.




Western blot analysis of RAP1GAP in rat Brain, HeLa lysates using RAP1GAP antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using RAP1GAP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: RAP1GAP Rabbit Monoclonal Antibody
产品货号: AMRe02522



Immunocytochemistry analysis of RAP1GAP(green) in HeLa using RAP1GAP antibody, and DAPI(blue).

注意事项 (Note)

For research use only .