

**产品名称: RAP1GAP Rabbit Monoclonal Antibody**  
**产品货号: AMRe02523**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | RAP1GAP Rabbit Monoclonal Antibody     |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB, ICC/IF, IP                         |
| 种属反应性 (Reactivity)     | Human, Mouse, Rat                      |

## 产品性能 (Performance)

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                         |
| 修饰 (Modification)   | Unmodified                                                                                           |
| 同种型 (Isotype)       | IgG                                                                                                  |
| 克隆 (Clonality)      | Monoclonal                                                                                           |
| 形式 (Form)           | Liquid                                                                                               |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| 储存溶液 (Buffer)       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| 纯化方式 (Purification) | Affinity Purification                                                                                |

## 免疫原信息 (Immunogen)

|                        |                                                                                 |
|------------------------|---------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | RAP1GAP                                                                         |
| 别名 (Alternative Names) | RAP1GAP; KIAA0474; RAP1GA1; Rap1 GTPase-activating protein 1; Rap1GAP; Rap1GAP1 |
| 基因 ID (Gene ID)        | 5909                                                                            |
| 蛋白 ID (SwissProt ID)   | P47736.                                                                         |

## 产品应用 (Application)

|                          |                                                  |
|--------------------------|--------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 73 kDa; Observed MW: 95 kDa       |

## 研究背景 (Background)

产品名称: RAP1GAP Rabbit Monoclonal Antibody  
产品货号: AMRe02523




GTPase activator for the nuclear Ras-related regulatory protein RAP-1A (KREV-1), converting it to the putatively inactive GDP-bound state.


## 研究领域 (Research Area)

Signal Transduction


## 图片 (Image Data)



Western blot analysis of RAP1GAP in A549 lysates using RAP1GAP antibody.



Western blot analysis of RAP1GAP in C6, 3T3 lysates using RAP1GAP antibody.



Immunocytochemistry analysis of RAP1GAP(green) in HeLa using RAP1GAP antibody, and DAPI(blue)

产品名称: RAP1GAP Rabbit Monoclonal Antibody  
产品货号: AMRe02523

---



**注意事项 (Note)**

For research use only .