

产品名称: GTPase HRAS Rabbit Monoclonal Antibody
产品货号: AMRe02525



产品概述 (Summary)

产品名称 (Production Name)	GTPase HRAS Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NRAS
别名 (Alternative Names)	NS6; CMNS; NCMS; ALPS4; N-ras; NRAS1
基因 ID (Gene ID)	4893
蛋白 ID (SwissProt ID)	P01111.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 21 kDa

研究背景 (Background)

Ras proteins bind GDP/GTP and possess intrinsic GTPase activity. Miscellaneous Mutations which change AA 12, 13 or 61

产品名称: GTPase HRAS Rabbit Monoclonal Antibody
产品货号: AMRe02525




activate the potential of Ras to transform cultured cells and are implicated in a variety of human tumors.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of GTPase HRAS in mouse brain lysates using GTPase HRAS antibody.



Western blot analysis of RAS in rat Brain lysates using RAS antibody.

注意事项 (Note)

For research use only .