

产品名称: RASSF2 Rabbit Monoclonal Antibody
产品货号: AMRe02526



产品概述 (Summary)

产品名称 (Production Name)	RASSF2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RASSF2
别名 (Alternative Names)	CENP-34; RASFADIN
基因 ID (Gene ID)	9770
蛋白 ID (SwissProt ID)	P50749.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa; Observed MW: 38 kDa

研究背景 (Background)

Potential tumor suppressor. Acts as a KRAS-specific effector protein. May promote apoptosis and cell cycle arrest. Stabilizes

产品名称: RASSF2 Rabbit Monoclonal Antibody
产品货号: AMRe02526




STK3/MST2 by protecting it from proteasomal degradation.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of RASSF2 in HeLa lysates using RASSF2 antibody.

注意事项 (Note)

For research use only .