

产品名称: RENT1 Rabbit Monoclonal Antibody
产品货号: AMRe02537



产品概述 (Summary)

产品名称 (Production Name)	RENT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UPF1 UPF1; KIAA0221; RENT1; Regulator of nonsense transcripts 1; ATP-dependent
别名 (Alternative Names)	helicase RENT1; Nonsense mRNA reducing factor 1; NORF1; Up-frameshift suppressor 1 homolog; hUpf1
基因 ID (Gene ID)	5976
蛋白 ID (SwissProt ID)	Q92900.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 124 kDa; Observed MW: 124 kDa

产品名称: RENT1 Rabbit Monoclonal Antibody
产品货号: AMRe02537




研究背景 (Background)

Plays a role in replication-dependent histone mRNA degradation at the end of phase S. Part of a post-splicing multiprotein complex. Involved in nonsense-mediated decay (NMD) as part of the SMG1C complex, a mRNA surveillance complex that recognizes and degrades mRNAs containing premature translation termination codons (PTCs).


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)




Western blot analysis of RENT1 / hUPF1 in Jurkat, HeLa lysates using RENT1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using RENT1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: RENT1 Rabbit Monoclonal Antibody
产品货号: AMRe02537



Immunocytochemistry analysis of RENT1 (green) in HeLa using RENT1 antibody, and DAPI (blue)

注意事项 (Note)

For research use only .