

产品名称: Retinoid X Receptor alpha Rabbit Monoclonal Antibody
产品货号: AMRe02539



产品概述 (Summary)

产品名称 (Production Name)	Retinoid X Receptor alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RXRA
别名 (Alternative Names)	NR2B1; RXR alpha1; Rxra; RXRalpha1
基因 ID (Gene ID)	6256
蛋白 ID (SwissProt ID)	P19793.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 51 kDa; Observed MW: 51 kDa

研究背景 (Background)


产品名称: Retinoid X Receptor alpha Rabbit Monoclonal Antibody
产品货号: AMRe02539

Receptor for retinoic acid. Retinoic acid receptors bind as heterodimers to their target response elements in response to their ligands, all-trans or 9-cis retinoic acid, and regulate gene expression in various biological processes.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Retinoid X Receptor alpha in HeLa, C6 lysates using Retinoid X Receptor alpha antibody.



Immunocytochemistry analysis of Retinoid X Receptor alpha (green) in MCF-7 using Retinoid X Receptor alpha antibody, and DAPI(blue)

注意事项 (Note)

For research use only .