

产品名称: Complex III Subunit 5 Rabbit Monoclonal Antibody
产品货号: AMRe02543



产品概述 (Summary)

产品名称 (Production Name)	Complex III Subunit 5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UQCRC1
别名 (Alternative Names)	petC; PGR1; RIP1; RIS1; RISP; UQCRC5; UQCRC1
基因 ID (Gene ID)	7386
蛋白 ID (SwissProt ID)	P47985.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 30 kDa; Observed MW: 24 kDa

研究背景 (Background)

产品名称: Complex III Subunit 5 Rabbit Monoclonal Antibody
产品货号: AMRe02543




Component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is a respiratory chain that generates an electrochemical potential coupled to ATP synthesis.


研究领域 (Research Area)

Microbiology


图片 (Image Data)



Western blot analysis of RISP in K562, rat Brain, 3T3, HeLa lysates using Complex III Subunit 5 antibody.



Immunocytochemistry analysis of RISP(green) in HeLa using RISP antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using UQCRCF51/RISP antibody. High-pressure

产品名称: Complex III Subunit 5 Rabbit Monoclonal Antibody
产品货号: AMRe02543



and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .