

产品名称: RBX1 Rabbit Monoclonal Antibody
产品货号: AMRe02546



产品概述 (Summary)

产品名称 (Production Name)	RBX1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBX1
别名 (Alternative Names)	ROC1; RNF75; BA554C12.1
基因 ID (Gene ID)	9978
蛋白 ID (SwissProt ID)	P62877.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 12 kDa; Observed MW: 12 kDa

研究背景 (Background)

E3 ubiquitin ligase component of multiple cullin-RING-based E3 ubiquitin-protein ligase complexes which mediate the

产品名称: RBX1 Rabbit Monoclonal Antibody
产品货号: AMRe02546




ubiquitination and subsequent proteasomal degradation of target proteins, including proteins involved in cell cycle progression, signal transduction, transcription and transcription-coupled nucleotide excision repair.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of ROC1 in Jurkat, THP-1, C6, 3T3, HeLa lysates using RBX1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using ROC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of ROC1 (green) in HeLa using ROC1 antibody, and DAPI (blue)

产品名称: RBX1 Rabbit Monoclonal Antibody
产品货号: AMRe02546



注意事项 (Note)

For research use only .