

产品名称: RPS6 Rabbit Monoclonal Antibody
产品货号: AMRe02553



产品概述 (Summary)

产品名称 (Production Name)	RPS6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPS6
别名 (Alternative Names)	RPS6; OK/SW-cl.2; 40S ribosomal protein S6; Phosphoprotein NP33
基因 ID (Gene ID)	6194
蛋白 ID (SwissProt ID)	P62753.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 29 kDa

研究背景 (Background)

May play an important role in controlling cell growth and proliferation through the selective translation of particular classes

产品名称: RPS6 Rabbit Monoclonal Antibody
产品货号: AMRe02553




of mRNA.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using RPS6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Western blot analysis of S6 Ribosomal Protein in Ramos, Jurkat, C6, 3T3 and HeLa lysates using S6 Ribosomal Protein antibody.




Immunocytochemistry analysis of HeLa using anti-S6 Ribosomal Protein antibody.

产品名称: RPS6 Rabbit Monoclonal Antibody
产品货号: AMRe02553



Immunoprecipitation analysis of HeLa cell lysates using S6 Ribosomal Protein antibody.



Western blot analysis of S6 Ribosomal Protein in CHO-K1, Raji, Raw264.7 and COS7 lysates using S6 Ribosomal Protein antibody.

注意事项 (Note)

For research use only .