

产品名称: RRM2 Rabbit Monoclonal Antibody
产品货号: AMRe02556



产品概述 (Summary)

产品名称 (Production Name)	RRM2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RRM2
别名 (Alternative Names)	R2; RR2M
基因 ID (Gene ID)	6241
蛋白 ID (SwissProt ID)	P31350.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 45 kDa

研究背景 (Background)

Provides the precursors necessary for DNA synthesis. Catalyzes the biosynthesis of deoxyribonucleotides from the

产品名称: RRM2 Rabbit Monoclonal Antibody
产品货号: AMRe02556




corresponding ribonucleotides. Inhibits Wnt signaling.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Immunocytochemistry analysis of RRM2 (green) in Jurkat using RRM2 antibody, and DAPI (blue).



Western blot analysis of RRM2 in K562, HeLa lysates using RRM2 antibody.

注意事项 (Note)

For research use only .