

产品名称: RUNX Rabbit Monoclonal Antibody
产品货号: AMRe02562



产品概述 (Summary)

产品名称 (Production Name)	RUNX Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RUNX1/RUNX2/RUNX3
别名 (Alternative Names)	RUNX1; AML1; RUNX2; OSF2; RUNX3; CBFA3
基因 ID (Gene ID)	861/864/860
蛋白 ID (SwissProt ID)	Q01196/Q13761/Q13950.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 49 kDa; Observed MW: 49 kDa

研究背景 (Background)

RUNX2 regulates the transcription of various genes including osteopontin, bone sialoprotein, and osteocalcin via binding

产品名称: RUNX Rabbit Monoclonal Antibody
产品货号: AMRe02562




to the core site of the enhancers or promoters. RUNX3/AML2 is a member of the Runt family of transcription factors. RUNX3 is important for the suppression of cell proliferation in the gastric epithelium, neurogenesis of the dorsal root ganglia, and T cell differentiation.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded mouse colon using RUNX antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of RUNX1/2/3 in RAW264.7 lysates using RUNX1/2/3 antibody.

注意事项 (Note)

For research use only .