

产品名称: Sall4 Rabbit Monoclonal Antibody
产品货号: AMRe02570



产品概述 (Summary)

产品名称 (Production Name)	Sall4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SALL4
别名 (Alternative Names)	Zinc finger protein 797; Zinc finger protein SALL4
基因 ID (Gene ID)	57167
蛋白 ID (SwissProt ID)	Q9UJQ4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 112 kDa; Observed MW: 112 kDa

研究背景 (Background)

Transcription factor with a key role in the maintenance and self-renewal of embryonic and hematopoietic stem cells.


产品名称: Sall4 Rabbit Monoclonal Antibody
产品货号: AMRe02570



研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Sall4 in Daudi lysates using Sall4 antibody.

注意事项 (Note)

For research use only .