

产品名称: Sarcomeric alpha Actinin Rabbit Monoclonal Antibody
产品货号: AMRe02572



产品概述 (Summary)

产品名称 (Production Name)	Sarcomeric alpha Actinin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACTN2
别名 (Alternative Names)	Actin binding protein; Actinin alpha 2; ACTN2; Alpha actinin 2; CMD1AA
基因 ID (Gene ID)	88
蛋白 ID (SwissProt ID)	P35609.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 104 kDa; Observed MW: 104 kDa

研究背景 (Background)

产品名称: Sarcomeric alpha Actinin Rabbit Monoclonal Antibody
产品货号: AMRe02572




F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded mouse colon using Sarcomeric alpha Actinin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Sarcomeric Alpha Actinin in mouse skeletal muscle, mouse heart lysates using Sarcomeric Alpha Actinin antibody.

产品名称: Sarcomeric alpha Actinin Rabbit Monoclonal Antibody
产品货号: AMRe02572



Western blot analysis of Sarcomeric Alpha Actinin in rat Brain lysates using Sarcomeric Alpha Actinin antibody

注意事项 (Note)

For research use only .