

产品名称: SCAI Rabbit Monoclonal Antibody
产品货号: AMRe02574



产品概述 (Summary)

产品名称 (Production Name)	SCAI Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SCAI
别名 (Alternative Names)	NET40; C9orf126
基因 ID (Gene ID)	286205
蛋白 ID (SwissProt ID)	Q8N9R8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 60 kDa

研究背景 (Background)

Tumor suppressor which functions to suppress MRTFA-induced SRF transcriptional activity. May function in the RHOA-

产品名称: SCAI Rabbit Monoclonal Antibody
产品货号: AMRe02574




DIAPH1 signal transduction pathway and regulate cell migration through transcriptional regulation of ITGB1.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of SCAI in C6, 3T3 lysates using SCAI antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SCAI antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .