

产品名称: SET Rabbit Monoclonal Antibody
产品货号: AMRe02590



产品概述 (Summary)

产品名称 (Production Name)	SET Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SET
别名 (Alternative Names)	2PP2A; IGAAD; TAF-I; I2PP2A; IPP2A2; PHAPII; TAF-IBETA
基因 ID (Gene ID)	6418
蛋白 ID (SwissProt ID)	Q01105.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 39 kDa

研究背景 (Background)

Multitasking protein, involved in apoptosis, transcription, nucleosome assembly and histone binding. Isoform 2 anti-

产品名称: SET Rabbit Monoclonal Antibody
产品货号: AMRe02590




apoptotic activity is mediated by inhibition of the GZMA-activated DNase, NME1.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Western blot analysis of SET in K562, HeLa lysates using SET antibody.



Immunocytochemistry analysis of SET(green) in HeLa using SET antibody,and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SET antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: SET Rabbit Monoclonal Antibody
产品货号: AMRe02590



注意事项 (Note)

For research use only .