

产品名称: SF3B3 Rabbit Monoclonal Antibody
产品货号: AMRe02592



产品概述 (Summary)

产品名称 (Production Name)	SF3B3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SF3B3 SF3B3; KIAA0017; SAP130; Splicing factor 3B subunit 3; Pre-mRNA-splicing
别名 (Alternative Names)	factor SF3b 130 kDa subunit; SF3b130; STAF130; Spliceosome-associated protein 130; SAP 130
基因 ID (Gene ID)	23450
蛋白 ID (SwissProt ID)	Q15393.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 136 kDa; Observed MW: 136 kDa

产品名称: SF3B3 Rabbit Monoclonal Antibody
产品货号: AMRe02592




研究背景 (Background)

Subunit of the splicing factor SF3B required for 'A' complex assembly formed by the stable binding of U2 snRNP to the branchpoint sequence (BPS) in pre-mRNA. Sequence independent binding of SF3A/SF3B complex upstream of the branch site is essential, it may anchor U2 snRNP to the pre-mRNA.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of SF3B3 in Raw264.7 lysates using SF3B3 antibody.



Western blot analysis of SF3B3 in SH-SY5Y, C6 lysates using SF3B3 antibody.



产品名称: SF3B3 Rabbit Monoclonal Antibody
产品货号: AMRe02592



Immunocytochemistry analysis of SF3B3(green) in Hela using SF3B3 antibody,and DAPI(blue)

注意事项 (Note)

For research use only .