

产品名称: SFRP4 Rabbit Monoclonal Antibody
产品货号: AMRe02594



产品概述 (Summary)

产品名称 (Production Name)	SFRP4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SFRP4
别名 (Alternative Names)	FRP4; FrpHE; PYL; Sfrp4
基因 ID (Gene ID)	6424
蛋白 ID (SwissProt ID)	Q6FHJ7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa; Observed MW: 53-55 kDa

研究背景 (Background)

Soluble frizzled-related proteins (sFRPS) function as modulators of Wnt signaling through direct interaction with Wnts.

产品名称: SFRP4 Rabbit Monoclonal Antibody
产品货号: AMRe02594




They have a role in regulating cell growth and differentiation in specific cell types.


研究领域 (Research Area)

Stem Cells

图片 (Image Data)



Western blot analysis of SFRP4 in rat kidney, rat heart lysates using SFRP4 antibody.



Western blot analysis of SFRP4 in K562, 3T3, HeLa lysates using SFRP4 antibody.

注意事项 (Note)

For research use only .