

产品名称: SNX9 Rabbit Monoclonal Antibody
产品货号: AMRe02597



产品概述 (Summary)

产品名称 (Production Name)	SNX9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNX9
别名 (Alternative Names)	SDP1; WISP; SH3PX1; SH3PXD3A
基因 ID (Gene ID)	51429
蛋白 ID (SwissProt ID)	Q9Y5X1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 67 kDa; Observed MW: 67 kDa

研究背景 (Background)

Involved in endocytosis and intracellular vesicle trafficking, both during interphase and at the end of mitosis.


产品名称: SNX9 Rabbit Monoclonal Antibody
产品货号: AMRe02597




研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of SNX9 in mouse intestines lysates using SNX9 antibody.




Western blot analysis of SH3PX1 in K562, Hela lysates using SH3PX1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using SNX9 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: SNX9 Rabbit Monoclonal Antibody
产品货号: AMRe02597



Immunocytochemistry analysis of SNX9 (green) in HeLa using SNX9 antibody, and DAPI (blue).

注意事项 (Note)

For research use only .