

产品名称: SMARCA5 Rabbit Monoclonal Antibody
产品货号: AMRe02620



产品概述 (Summary)

产品名称 (Production Name)	SMARCA5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMARCA5
别名 (Alternative Names)	ISWI; SNF2H; hISWI; hSNF2H; WCRF135
基因 ID (Gene ID)	8467
蛋白 ID (SwissProt ID)	O60264.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 122 kDa; Observed MW: 122 kDa

研究背景 (Background)

Helicase that possesses intrinsic ATP-dependent nucleosome-remodeling activity. Complexes containing SMARCA5 are

产品名称: SMARCA5 Rabbit Monoclonal Antibody
产品货号: AMRe02620




capable of forming ordered nucleosome arrays on chromatin; this may require intact histone H4 tails.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Immunocytochemistry analysis of SMARCA5 (green) in Jurkat using SMARCA5 antibody, and DAPI (blue).



Western blot analysis of SNF2H in K562, C6, 3T3, HeLa lysates using SNF2H antibody.

注意事项 (Note)

For research use only .