

产品名称: SMARCB1 Rabbit Monoclonal Antibody
产品货号: AMRe02622



产品概述 (Summary)

产品名称 (Production Name)	SMARCB1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMARCB1
别名 (Alternative Names)	SMARCB1; BAF47; hSNF5; INI1; RDT; RTPS1; Sfh1p; SMARCB1; SNF5 homolog; SNF5L1; Snr1; SWI/SNF comp
基因 ID (Gene ID)	6598
蛋白 ID (SwissProt ID)	Q12824.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa; Observed MW: 44 kDa

研究背景 (Background)

产品名称: SMARCB1 Rabbit Monoclonal Antibody
产品货号: AMRe02622




The SWI-SNF complex is involved in the activation of transcription via the remodeling of nucleosome structure in an ATP-dependent manner. Brm (also designated SNF2 α) and Brg-1 (also designated SNF2 β) are the ATPase subunits of the mammalian SWI-SNF complex.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using SMARCB1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of SNF5 (green) in HeLa using SNF5 antibody, and DAPI (blue)

注意事项 (Note)

For research use only .