

产品名称: Spermine Synthase Rabbit Monoclonal Antibody
产品货号: AMRe02630



产品概述 (Summary)

产品名称 (Production Name)	Spermine Synthase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMS
别名 (Alternative Names)	SRS; SpS; MRSR; SPMSY
基因 ID (Gene ID)	6611
蛋白 ID (SwissProt ID)	P52788.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 41 kDa; Observed MW: 41 kDa

研究背景 (Background)

This gene encodes a protein belonging to the spermidine/spermin synthase family. Pseudogenes of this gene are located

产品名称: Spermine Synthase Rabbit Monoclonal Antibody
产品货号: AMRe02630




on chromosomes 1, 5, 6 and X. Mutations in this gene are associated with X-linked Snyder-Robinson mental retardation syndrome. Multiple transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Spermine synthase in K562, HeLa lysates using Spermine Synthase antibody.

注意事项 (Note)

For research use only .