

产品名称: SPINK5 Rabbit Monoclonal Antibody
产品货号: AMRe02631



产品概述 (Summary)

产品名称 (Production Name)	SPINK5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SPINK5
别名 (Alternative Names)	NS; NETS; LEKTI; LETKI; VAKTI
基因 ID (Gene ID)	11005
蛋白 ID (SwissProt ID)	Q9NQ38.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 121 kDa; Observed MW: 121 kDa

研究背景 (Background)

Serine protease inhibitor, probably important for the anti-inflammatory and/or antimicrobial protection of mucous

产品名称: SPINK5 Rabbit Monoclonal Antibody
产品货号: AMRe02631




epithelia. Contribute to the integrity and protective barrier function of the skin by regulating the activity of defense-activating and desquamation-involved proteases. Inhibits KLK5, it's major target, in a pH-dependent manner. Inhibits KLK7, KLK14 CASP14, and trypsin.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of SPINK5/LEKTI in HeLa lysates using SPINK5 antibody.

注意事项 (Note)

For research use only .