

产品名称: SRP72 Rabbit Monoclonal Antibody
产品货号: AMRe02635



产品概述 (Summary)

产品名称 (Production Name)	SRP72 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SRP72
别名 (Alternative Names)	BMFF; BMFS1; HEL103
基因 ID (Gene ID)	6731
蛋白 ID (SwissProt ID)	O76094.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 75 kDa; Observed MW: 75 kDa

研究背景 (Background)

Signal-recognition-particle assembly has a crucial role in targeting secretory proteins to the rough endoplasmic reticulum

产品名称: SRP72 Rabbit Monoclonal Antibody
产品货号: AMRe02635




membrane. Binds the 7S RNA only in presence of SRP68. This ribonucleoprotein complex might interact directly with the docking protein in the ER membrane and possibly participate in the elongation arrest function.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of SRP72 in HeLa, HUVEC, K562 lysates using SRP72 antibody.

注意事项 (Note)

For research use only .