

产品名称: Striatin 4 Rabbit Monoclonal Antibody
产品货号: AMRe02645



产品概述 (Summary)

产品名称 (Production Name)	Striatin 4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STRN4
别名 (Alternative Names)	ZIN; PPP2R6C; zinedin
基因 ID (Gene ID)	29888
蛋白 ID (SwissProt ID)	Q9NRL3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 81 kDa; Observed MW: 81 kDa

研究背景 (Background)

Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein. MiscellaneousThe name

产品名称: Striatin 4 Rabbit Monoclonal Antibody
产品货号: AMRe02645




'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane).

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Striatin 4 in HeLa lysates using Striatin 4 antibody.

注意事项 (Note)

For research use only .