

产品名称: STUB1 Rabbit Monoclonal Antibody
产品货号: AMRe02648



产品概述 (Summary)

产品名称 (Production Name)	STUB1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STUB1
别名 (Alternative Names)	STUB1; CHIP; HSPABP2; NY-CO-7; SDCCAG7; UBOX1
基因 ID (Gene ID)	10273
蛋白 ID (SwissProt ID)	Q9UNE7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 35 kDa; Observed MW: 35 kDa

研究背景 (Background)

E3 ubiquitin-protein ligase which targets misfolded chaperone substrates towards proteasomal degradation. Ubiquitinates

产品名称: STUB1 Rabbit Monoclonal Antibody
产品货号: AMRe02648




NOS1 in concert with Hsp70 and Hsp40. Modulates the activity of several chaperone complexes, including Hsp70, Hsc70 and Hsp90.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of STUB1 in K562, rat Brain, C6, 3T3, HeLa lysates using STUB1 antibody.



Immunocytochemistry analysis of STUB1(green) in HeLa using STUB1 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .