

产品名称: TCEA1 Rabbit Monoclonal Antibody
产品货号: AMRe02674



产品概述 (Summary)

产品名称 (Production Name)	TCEA1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TCEA1
别名 (Alternative Names)	SII; TCEA; TF2S; GTF2S; TFIIIS
基因 ID (Gene ID)	6917
蛋白 ID (SwissProt ID)	P23193.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 34 kDa

研究背景 (Background)

Necessary for efficient RNA polymerase II transcription elongation past template-encoded arresting sites. The arresting

产品名称: TCEA1 Rabbit Monoclonal Antibody
产品货号: AMRe02674




sites in DNA have the property of trapping a certain fraction of elongating RNA polymerases that pass through, resulting in locked ternary complexes. Cleavage of the nascent transcript by S-II allows the resumption of elongation from the new 3'-terminus. Miscellaneous S-II binds to RNA-polymerase II in the absence of transcription.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of TCEA1 in rat spleen lysates using TCEA1 antibody.



Western blot analysis of TCEA1 in EAHY, K562 lysates using TCEA1 antibody.

注意事项 (Note)

For research use only .