

产品名称: Transcriptional Enhancer Factor TEF 1 Rabbit
Monoclonal Antibody
产品货号: AMRe02680



产品概述 (Summary)

产品名称 (Production Name)	Transcriptional Enhancer Factor TEF 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TEAD1
别名 (Alternative Names)	TEAD1; TCF13; TEF1; Transcriptional enhancer factor TEF-1; NTEF-1; Protein GT-IIC; TEA domain family member 1; TEAD-1; Transcription factor 13; TCF-13
基因 ID (Gene ID)	7003
蛋白 ID (SwissProt ID)	P28347.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 48 kDa

产品名称: Transcriptional Enhancer Factor TEF 1 Rabbit
Monoclonal Antibody
产品货号: AMRe02680




研究背景 (Background)

Transcription factor which plays a key role in the Hippo signaling pathway, a pathway involved in organ size control and tumor suppression by restricting proliferation and promoting apoptosis.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of TEF1/TEAD1 in C6, CHO-K1, HeLa lysates using Transcriptional Enhancer Factor TEF 1 antibody.

注意事项 (Note)

For research use only .