

产品名称: Transmembrane Protein 9B Rabbit Monoclonal Antibody
产品货号: AMRe02700



产品概述 (Summary)

产品名称 (Production Name)	Transmembrane Protein 9B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TMEM9B
别名 (Alternative Names)	C11orf15
基因 ID (Gene ID)	56674
蛋白 ID (SwissProt ID)	Q9NQ34.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa

研究背景 (Background)

产品名称: Transmembrane Protein 9B Rabbit Monoclonal Antibody
产品货号: AMRe02700




Enhances production of proinflammatory cytokines induced by TNF, IL1B, and TLR ligands. Has a role in TNF activation of both the NF-kappaB and MAPK pathways.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of TMEM9B in Hacat, PC-12 lysates using Transmembrane Protein 9B antibody.



Western blot analysis of TMEM9B in 3T3 lysates using TMEM9B antibody.

注意事项 (Note)

For research use only .