

产品名称: TNFAIP3 Rabbit Monoclonal Antibody
产品货号: AMRe02701



产品概述 (Summary)

产品名称 (Production Name)	TNFAIP3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNFAIP3 TNFAIP3; OTUD7C; Tumor necrosis factor alpha-induced protein 3; TNF alpha-induced protein 3; OTU domain-containing protein 7C; Putative DNA-binding protein A20; Zinc finger protein A20
别名 (Alternative Names)	
基因 ID (Gene ID)	7128
蛋白 ID (SwissProt ID)	P21580.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 90 kDa; Observed MW: 82 kDa

产品名称: TNFAIP3 Rabbit Monoclonal Antibody
产品货号: AMRe02701



研究背景 (Background)

TNFAIP3 is cytokine-inducible protein that functions to inhibit apoptosis and activate NF- κ B. Play a role in the function of the lymphoid system. Required for LPS-induced production of proinflammatory cytokines and IFN β in LPS-tolerized macrophages.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Immunocytochemistry analysis of TNFAIP3 (green) in HL-60 using TNFAIP3 antibody, and DAPI (blue).



Western blot analysis of TNFAIP3 in HeLa lysates using TNFAIP3 antibody.



产品名称: TNFAIP3 Rabbit Monoclonal Antibody
产品货号: AMRe02701



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using TNFAIP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .