

产品名称: Tensin 4 Rabbit Monoclonal Antibody
产品货号: AMRe02702



产品概述 (Summary)

产品名称 (Production Name)	Tensin 4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNS4
别名 (Alternative Names)	CTEN; C-terminal tensin-like protein
基因 ID (Gene ID)	84951
蛋白 ID (SwissProt ID)	Q8IZW8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 77 kDa; Observed MW: 90 kDa

研究背景 (Background)

May be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton .

产品名称: Tensin 4 Rabbit Monoclonal Antibody
产品货号: AMRe02702




May promote apoptosis, via its cleavage by caspase-3.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of TNS4/CTEN in Raw264.7, C6 lysates using Tensin 4 antibody.



Western blot analysis of TNS4/CTEN in K562 lysates using TNS4/CTEN antibody.

注意事项 (Note)

For research use only .