

产品名称: Transferrin Receptor 1 Rabbit Monoclonal Antibody
产品货号: AMRe02710



产品概述 (Summary)

产品名称 (Production Name)	Transferrin Receptor 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TFRC
别名 (Alternative Names)	TFRC; Transferrin receptor protein 1; TR; TfR; TfR1; Trfr; T9; p90; CD71
基因 ID (Gene ID)	7037
蛋白 ID (SwissProt ID)	P02786.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 85 kDa; Observed MW: 90 kDa

研究背景 (Background)

产品名称: Transferrin Receptor 1 Rabbit Monoclonal Antibody
产品货号: AMRe02710




Binds the telomeric double-stranded TTAGGG repeat and negatively regulates telomere length. Involved in the regulation of the mitotic spindle. Component of the shelterin complex (telosome) that is involved in the regulation of telomere length and protection.


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Transferrin Receptor in K562, HeLa lysates using Transferrin Receptor 1 antibody.



Immunocytochemistry analysis of Transferrin Receptor(green) in HeLa using Transferrin Receptor antibody, and DAPI(blue)

注意事项 (Note)

For research use only .