

产品名称: Translin Associated Protein X Rabbit Monoclonal Antibody
产品货号: AMRe02712



产品概述 (Summary)

产品名称 (Production Name)	Translin Associated Protein X Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TSNAX
别名 (Alternative Names)	TRAX
基因 ID (Gene ID)	7257
蛋白 ID (SwissProt ID)	Q99598.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 33 kDa

研究背景 (Background)

产品名称: Translin Associated Protein X Rabbit Monoclonal Antibody
产品货号: AMRe02712




Acts in combination with TSN as an endonuclease involved in the activation of the RNA-induced silencing complex (RISC).
Possible role in spermatogenesis.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of TRAX in 3T3, mouse testis lysates using Translin Associated Protein X antibody.



Western blot analysis of TRAX in K562, rat Brain, C6, HeLa lysates using TRAX antibody.

注意事项 (Note)

For research use only .