

产品名称: UCHL5 Rabbit Monoclonal Antibody
产品货号: AMRe02740



产品概述 (Summary)

产品名称 (Production Name)	UCHL5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UCHL5
别名 (Alternative Names)	UCH37; CGI-70; INO80R; UCH-L5
基因 ID (Gene ID)	51377
蛋白 ID (SwissProt ID)	Q9Y5K5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa; Observed MW: 38 kDa

研究背景 (Background)

Protease that specifically cleaves 'Lys-48'-linked polyubiquitin chains.


产品名称: UCHL5 Rabbit Monoclonal Antibody
产品货号: AMRe02740




研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of UCH37 in rat spleen lysates using UCHL5 antibody.



Western blot analysis of UCH37 in 3T3, Hela lysates using UCH37 antibody.



Immunocytochemistry analysis of UCH37(green) in Hela using UCH37 antibody, and DAPI(blue)

注意事项 (Note)

产品名称: UCHL5 Rabbit Monoclonal Antibody
产品货号: AMRe02740



For research use only .