

产品名称: UFC1 Rabbit Monoclonal Antibody
产品货号: AMRe02743



产品概述 (Summary)

产品名称 (Production Name)	UFC1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UFC1
别名 (Alternative Names)	NEDSG; HSPC155
基因 ID (Gene ID)	51506
蛋白 ID (SwissProt ID)	Q9Y3C8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 19 kDa; Observed MW: 19 kDa

研究背景 (Background)

E2-like enzyme which forms an intermediate with UFM1 via a thioester linkage.


产品名称: UFC1 Rabbit Monoclonal Antibody
产品货号: AMRe02743




研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of UFC1 in K562, HeLa lysates using UFC1 antibody.



Immunocytochemistry analysis of UFC1(green) in HeLa using UFC1 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .