

产品名称: Vitronectin Rabbit Monoclonal Antibody
产品货号: AMRe02765



产品概述 (Summary)

产品名称 (Production Name)	Vitronectin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VTN
别名 (Alternative Names)	VTN; Vitronectin; VN; S-protein; Serum-spreading factor; V75
基因 ID (Gene ID)	7448
蛋白 ID (SwissProt ID)	P04004.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 70 kDa

研究背景 (Background)

The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell

产品名称: Vitronectin Rabbit Monoclonal Antibody
产品货号: AMRe02765




adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Vitronectin in Human serum lysates using Vitronectin antibody.

注意事项 (Note)

For research use only .