

**产品名称:** Phospho-c Myb (Ser11) Rabbit Monoclonal Antibody

**产品货号:** AMRe02766

### 产品概述 (Summary)

|                   |                                                                                                      |
|-------------------|------------------------------------------------------------------------------------------------------|
| 描述(Description)   | Recombinant rabbit monoclonal antibody                                                               |
| 宿主(Host)          | Rabbit                                                                                               |
| 应用(Application)   | WB,IHC,IP,ChIP                                                                                       |
| 种属反应(Reactivity)  | Human                                                                                                |
| 偶联物(Conjugation)  | Unconjugated                                                                                         |
| 修饰(Modification)  | Phosphorylated                                                                                       |
| 同种型(Isotype)      | IgG                                                                                                  |
| 克隆(Clonality)     | Monoclonal                                                                                           |
| 剂型(Form)          | Liquid                                                                                               |
| 保存条件(Storage)     | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| 储存溶液(Buffer)      | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| 纯化 (Purification) | Affinity Purification                                                                                |

### 产品应用(Application)

|                       |                                                       |
|-----------------------|-------------------------------------------------------|
| 稀释比(Dilution Ratio)   | WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50,ChIP 1:20 |
| 分子量(Molecular Weight) | Calculated MW: 72 kDa; Observed MW: 80 kDa            |

### 抗原信息(Antigen Information)

|                       |                                                                             |
|-----------------------|-----------------------------------------------------------------------------|
| 基因名(Gene Name)        | MYB                                                                         |
| 别名(Alternative Names) | MYB; Transcriptional activator Myb; Proto-oncogene c-Myb                    |
| 基因 ID(Gene ID)        | 4602                                                                        |
| SwissProt ID          | P10242                                                                      |
| 免疫原(Immunogen)        | A synthetic Phosphorylated peptide corresponding to residues target protein |

### 研究背景 (Background)

Transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC[GT]G-3'. Plays an important role in the control of proliferation and differentiation of hematopoietic progenitor cells.


## 研究领域 (Research Area)

Epigenetics and Nuclear Signaling


## 注意事项 (Note)

For research use only.

## 图片 (Image Data)



Western blot analysis of Phospho-vMyb + c-Myb (Ser11) in Jurkat lysates using Phospho-c Myb (Ser11) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using vMyb/c-Myb (Phospho-Ser11) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

## EnkiLife 优势产品推荐

| 产品类别                        | 产品货号     | 产品名称                                               |
|-----------------------------|----------|----------------------------------------------------|
| WB 解决方案                     | RA10020  | 2h 极速 WB 即用型全流程试剂盒                                 |
|                             | RA10021  | 4h 快速 WB 即用型全流程试剂盒                                 |
|                             | RA10042  | 5h 畅享版 WB 全流程试剂盒                                   |
|                             | RA10037  | 校准级彩色预染蛋白 Marker (8-180kDa)                        |
|                             | RA10038  | 校准级彩色预染蛋白 Marker (10-250kDa)                       |
|                             | RA10039  | 校准级高分子彩色预染蛋白 Marker (25-400kDa)                    |
| TSA 多重荧光染色试剂盒               | RA10008  | TSA 双标三色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10009  | TSA 三标四色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10010  | TSA 四标五色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10011  | TSA 五标六色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10012  | TSA 六标七色多重荧光染色试剂盒 (mIHC)                           |
| IHC 检测试剂盒                   | RA10006  | HRP Anti-Mouse/Rabbit IHC Detection System         |
|                             | RA10007  | Polymer-HRP Anti-Mouse/Rabbit IHC Detection System |
| 抗体标记试剂盒                     | RE80004p | 辣根过氧化物酶(HRP)抗体标记试剂盒                                |
|                             | RE80002q | Sulfo-NHS-生物素标记试剂盒                                 |
|                             | RE80007p | Cy3 荧光素标记试剂盒                                       |
|                             | RE80011p | Fluor488 荧光素标记试剂盒                                  |
|                             | RE80017p | Fluor750 荧光素标记试剂盒                                  |
|                             | RE80005p | 藻红蛋白(R-PE) 抗体快速标记试剂盒                               |
|                             | RE80040  | PE-Cy7 串联染料抗体快速标记试剂盒                               |
| 稳转细胞系构建服务<br>(免费赠送全膜 WB 验证) | TS-0001  | 过表达稳转细胞系构建                                         |
|                             | TS-0002  | 敲低稳转细胞系构建                                          |
|                             | TS-0003  | 敲除细胞系构建                                            |

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>