

产品名称: VPS41 Rabbit Monoclonal Antibody  
产品货号: AMRe02768



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | VPS41 Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB                                     |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                         |
| 修饰 (Modification)   | Unmodified                                                                                           |
| 同种型 (Isotype)       | IgG                                                                                                  |
| 克隆 (Clonality)      | Monoclonal                                                                                           |
| 形式 (Form)           | Liquid                                                                                               |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| 储存溶液 (Buffer)       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| 纯化方式 (Purification) | Affinity Purification                                                                                |

## 免疫原信息 (Immunogen)

|                        |                         |
|------------------------|-------------------------|
| 基因名 (Gene Name)        | VPS41                   |
| 别名 (Alternative Names) | HVPS41; HVSP41; hVps41p |
| 基因 ID (Gene ID)        | 27072                   |
| 蛋白 ID (SwissProt ID)   | P49754.                 |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000                            |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 99 kDa; Observed MW: 99 kDa |

## 研究背景 (Background)

Plays a role in vesicle-mediated protein trafficking to lysosomal compartments including the endocytic membrane

产品名称: VPS41 Rabbit Monoclonal Antibody  
产品货号: AMRe02768




transport and autophagic pathways. Believed to act in part as a core component of the putative HOPS endosomal tethering complex is proposed to be involved in the Rab5-to-Rab7 endosome conversion probably implicating MON1A/B, and via binding SNAREs and SNARE complexes to mediate tethering and docking events during SNARE-mediated membrane fusion. The HOPS complex is proposed to be recruited to Rab7 on the late endosomal membrane and to regulate late endocytic, phagocytic and autophagic traffic towards lysosomes (PubMed:23351085). Involved in homotypic vesicle fusions between late endosomes and in heterotypic fusions between late endosomes and lysosomes implicated in degradation of endocytosed cargo (PubMed:9159129, PubMed:23167963, PubMed:25445562, PubMed:25908847). Required for fusion of autophagosomes with lysosomes (PubMed:25783203). May link the HOPS complex to endosomal Rab7 via its association with RILP and to lysosomal membranes via its association with ARL8B, suggesting that these interactions may bring the compartments to close proximity for fusion (PubMed:25445562, PubMed:25908847). Involved in the direct trans-Golgi network to late endosomes transport of lysosomal membrane proteins independently of HOPS (PubMed:23322049). Involved in sorting to the regulated secretory pathway presumably implicating the AP-3 adaptor complex . May play a role in HOPS-independent function in the regulated secretory pathway (PubMed:24210660).


## 研究领域 (Research Area)

Signal Transduction

## 图片 (Image Data)



Western blot analysis of VPS41 in rat heart lysates using VPS41 antibody.



**产品名称: VPS41 Rabbit Monoclonal Antibody**  
**产品货号: AMRe02768**



---

Western blot analysis of VPS41 in Hela, 293T, mouse brain lysates using VPS41 antibody.

**注意事项 (Note)**

For research use only .