

产品名称: WAVE 1 Rabbit Monoclonal Antibody
产品货号: AMRe02770



产品概述 (Summary)

产品名称 (Production Name)	WAVE 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	WASF1
别名 (Alternative Names)	WAVE; SCAR1; WAVE1
基因 ID (Gene ID)	8936
蛋白 ID (SwissProt ID)	Q92558.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa; Observed MW: 90 kDa

研究背景 (Background)

Downstream effector molecule involved in the transmission of signals from tyrosine kinase receptors and small GTPases to

产品名称: WAVE 1 Rabbit Monoclonal Antibody
产品货号: AMRe02770




the actin cytoskeleton. Promotes formation of actin filaments. Part of the WAVE complex that regulates lamellipodia formation. The WAVE complex regulates actin filament reorganization via its interaction with the Arp2/3 complex . As component of the WAVE1 complex, required for BDNF-NTRK2 endocytic trafficking and signaling from early endosomes .

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of WAVE 1 in rat Brain, C6, 3T3, HeLa lysates using WAVE 1 antibody.

注意事项 (Note)

For research use only .