

产品名称: Wnt2b Rabbit Monoclonal Antibody
产品货号: AMRe02774



产品概述 (Summary)

产品名称 (Production Name)	Wnt2b Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	WNT2B
别名 (Alternative Names)	WNT13; Wnt2b; XWNT2
基因 ID (Gene ID)	7482
蛋白 ID (SwissProt ID)	Q93097.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa; Observed MW: 44 kDa

研究背景 (Background)

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a

产品名称: Wnt2b Rabbit Monoclonal Antibody
产品货号: AMRe02774




signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.


研究领域 (Research Area)

Stem Cells

图片 (Image Data)



Immunocytochemistry analysis of Wnt2b (green) in 293t using Wnt2b antibody, and DAPI (blue).



Western blot analysis of Wnt2b in K562, C6, 3T3, HeLa lysates using Wnt2b antibody.

注意事项 (Note)

For research use only .