

产品名称: XIAP Rabbit Monoclonal Antibody
产品货号: AMRe02777



产品概述 (Summary)

产品名称 (Production Name)	XIAP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Xiap
别名 (Alternative Names)	API3; ILP1; MIHA; XLP2; BIRC4; IAP-3; hIAP3; hIAP-3
基因 ID (Gene ID)	11798
蛋白 ID (SwissProt ID)	Q60989.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 56 kDa; Observed MW: 53 kDa

研究背景 (Background)

Multi-functional protein which regulates not only caspases and apoptosis, but also modulates inflammatory signaling and

产品名称: XIAP Rabbit Monoclonal Antibody
产品货号: AMRe02777




immunity, copper homeostasis, mitogenic kinase signaling, cell proliferation, as well as cell invasion and metastasis. Acts as a direct caspase inhibitor. Directly bind to the active site pocket of CASP3 and CASP7 and obstructs substrate entry. Inactivates CASP9 by keeping it in a monomeric, inactive state.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of XIAP in K562, rat Brain, 3T3 lysates using XIAP antibody.

注意事项 (Note)

For research use only .