

产品名称: YAP1 Rabbit Monoclonal Antibody
产品货号: AMRe02781



产品概述 (Summary)

产品名称 (Production Name)	YAP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YAP1
别名 (Alternative Names)	YAP1; YAP65; Yorkie homolog; 65 kDa Yes-associated protein; YAP65
基因 ID (Gene ID)	10413
蛋白 ID (SwissProt ID)	P46937.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 70-75 kDa

研究背景 (Background)

YAP (Yes-associated protein, YAP65) was identified based on its ability to associate with the SH3 domain of Yes.

产品名称: YAP1 Rabbit Monoclonal Antibody
产品货号: AMRe02781




Transcriptional regulator which can act both as a coactivator and a corepressor and is the critical downstream regulatory target in the Hippo signaling pathway that plays a pivotal role in organ size control and tumor suppression by restricting proliferation and promoting apoptosis.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Immunofluorescence analysis of YAP1 in mouse embryo using YAP1 antibody.



Western blot analysis of YAP1 in HeLa lysates using YAP1 antibody.

注意事项 (Note)

For research use only .