

产品名称: Tissue Factor Rabbit Monoclonal Antibody
产品货号: AMRe02797



产品概述 (Summary)

产品名称 (Production Name)	Tissue Factor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	F3
别名 (Alternative Names)	F3; Tissue factor; TF; Coagulation factor III; Thromboplastin; CD antigen CD142
基因 ID (Gene ID)	2152
蛋白 ID (SwissProt ID)	P13726.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 45 kDa

研究背景 (Background)

Initiates blood coagulation by forming a complex with circulating factor VII or VIIa. The [TF:VIIa] complex activates factors

产品名称: Tissue Factor Rabbit Monoclonal Antibody
产品货号: AMRe02797




IX or X by specific limited protolysis. TF plays a role in normal hemostasis by initiating the cell-surface assembly and propagation of the coagulation protease cascade.


研究领域 (Research Area)

Cardiovascular


图片 (Image Data)



Western blot analysis of Tissue Factor in mouse brain, rat brain lysates using Tissue Factor antibody.



Western blot analysis of Tissue Factor in Hela, A549, MCF-7, U251, HL-60 lysates using Tissue Factor antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Tissue Factor antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Tissue Factor Rabbit Monoclonal Antibody
产品货号: AMRe02797



注意事项 (Note)

For research use only .