

产品名称: TIM 1 Rabbit Monoclonal Antibody
产品货号: AMRe02798



产品概述 (Summary)

产品名称 (Production Name)	TIM 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HAVCR1
别名 (Alternative Names)	CD365; HAVCR 1; HAVCR; Havcr1; KIM 1; TIM; TIM1; TIMD1
基因 ID (Gene ID)	26762
蛋白 ID (SwissProt ID)	Q96D42.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa; Observed MW: 50 kDa

研究背景 (Background)

May play a role in T-helper cell development and the regulation of asthma and allergic diseases. Receptor for TIMD4 (By

产品名称: TIM 1 Rabbit Monoclonal Antibody
产品货号: AMRe02798




similarity). In case of human hepatitis A virus (HHAV) infection, functions as a cell-surface receptor for the virus. May play a role in kidney injury and repair.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of TIM 1 in Raw264.7, PC-12 lysates using Hepatitis A Virus ular Receptor 1 antibody.



Western blot analysis of TIM 1 in A549, MCF-7, U251, HL-60 lysates using TIM 1 antibody.

注意事项 (Note)

For research use only .