

产品名称: SF3A3 Rabbit Monoclonal Antibody
产品货号: AMRe02799



产品概述 (Summary)

产品名称 (Production Name)	SF3A3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SF3A3
别名 (Alternative Names)	splicing factor 3a subunit 3; PRP9; PRPF9; SAP61; SF3a60
基因 ID (Gene ID)	10946
蛋白 ID (SwissProt ID)	Q12874.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 59 kDa; Observed MW: 59 kDa

研究背景 (Background)

Involved in pre-mRNA splicing as a component of the splicing factor SF3A complex that contributes to the assembly of the

产品名称: SF3A3 Rabbit Monoclonal Antibody
产品货号: AMRe02799




17S U2 snRNP, and the subsequent assembly of the pre-spliceosome 'E' complex and the pre-catalytic spliceosome 'A' complex.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of SF3A3 in 3T3, Raw264.7, C6 lysates using SF3A3 antibody.



Western blot analysis of SF3A3 in HeLa, CHO-K1, C6, Jurkat lysates using SF3A3 antibody.

注意事项 (Note)

For research use only .