

产品名称: Hex Rabbit Monoclonal Antibody  
产品货号: AMRe02806



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | Hex Rabbit Monoclonal Antibody         |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IP                                  |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                         |
| 修饰 (Modification)   | Unmodified                                                                                           |
| 同种型 (Isotype)       | IgG                                                                                                  |
| 克隆 (Clonality)      | Monoclonal                                                                                           |
| 形式 (Form)           | Liquid                                                                                               |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| 储存溶液 (Buffer)       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| 纯化方式 (Purification) | Affinity Purification                                                                                |

## 免疫原信息 (Immunogen)

|                        |                            |
|------------------------|----------------------------|
| 基因名 (Gene Name)        | HHEX                       |
| 别名 (Alternative Names) | HEX; HHEX; HMPH; PRH; PRHX |
| 基因 ID (Gene ID)        | 3087                       |
| 蛋白 ID (SwissProt ID)   | Q03014.                    |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IP 1:20-1:50               |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 30 kDa; Observed MW: 30 kDa |

## 研究背景 (Background)

Recognizes the DNA sequence 5'-ATTAA-3'. Transcriptional repressor. May play a role in hematopoietic differentiation.

**产品名称: Hex Rabbit Monoclonal Antibody**  
**产品货号: AMRe02806**




Establishes anterior identity at two levels; acts early to enhance canonical WNT-signaling by repressing expression of TLE4, and acts later to inhibit NODAL-signaling by directly targeting NODAL.

### 研究领域 (Research Area)

Epigenetics and Nuclear Signaling

### 图片 (Image Data)



Western blot analysis of Hex in Jurkat, HeLa lysates using Hex antibody.

### 注意事项 (Note)

For research use only .