

产品名称: FIG4 Rabbit Monoclonal Antibody
产品货号: AMRe02815



产品概述 (Summary)

产品名称 (Production Name)	FIG4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FIG4
别名 (Alternative Names)	YVS; BTOP; SAC3; ALS11; CMT4J; KIAA0274; dJ249I4.1
基因 ID (Gene ID)	9896
蛋白 ID (SwissProt ID)	Q92562.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 104 kDa; Observed MW: 104 kDa

研究背景 (Background)

Dual specificity phosphatase component of the PI(3,5)P2 regulatory complex which regulates both the synthesis and


产品名称: FIG4 Rabbit Monoclonal Antibody
产品货号: AMRe02815



turnover of phosphatidylinositol 3,5-bisphosphate (PtdIns(3,5)P2).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FIG4 in HL-60, U2OS, U87-MG lysates using FIG4 antibody.

注意事项 (Note)

For research use only .