

产品名称: Histone H1.2 Rabbit Monoclonal Antibody
产品货号: AMRe02816



产品概述 (Summary)

产品名称 (Production Name)	Histone H1.2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	H1-2
别名 (Alternative Names)	H1.a; H1F2; H1s-1; HIST1H1C; Histone H1.2; Histone H1c; Histone H1d; Histone H1s-1
基因 ID (Gene ID)	3006
蛋白 ID (SwissProt ID)	P16403.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 30 kDa

研究背景 (Background)

产品名称: Histone H1.2 Rabbit Monoclonal Antibody
产品货号: AMRe02816




Histones H1 are necessary for the condensation of nucleosome chains into higher order structures.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Immunocytochemistry analysis of Histone H1.2 (green) in K562 using Histone H1.2 antibody, and DAPI (blue).



Western blot analysis of Histone H1.2 in 3T3 lysates using Histone H1.2 antibody.



Western blot analysis of Histone H1.2 in A549, MCF-7, U251, HL-60 lysates using Histone H1.2 antibody

产品名称: Histone H1.2 Rabbit Monoclonal Antibody
产品货号: AMRe02816



注意事项 (Note)

For research use only .