

产品名称: SF20 Rabbit Monoclonal Antibody
产品货号: AMRe02819



产品概述 (Summary)

产品名称 (Production Name)	SF20 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYDGF
别名 (Alternative Names)	IL25; IL27; SF20; IL27w; C19orf10; R33729_1; EUROIMAGE1875335
基因 ID (Gene ID)	56005
蛋白 ID (SwissProt ID)	Q969H8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 19 kDa; Observed MW: 19 kDa

研究背景 (Background)

The protein encoded by this gene was previously thought to support proliferation of lymphoid cells and was considered an

产品名称: SF20 Rabbit Monoclonal Antibody
产品货号: AMRe02819




interleukin. However, this activity has not been reproducible and the function of this protein is currently unknown.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of SF20 in HeLa, Jurkat lysates using SF20 antibody.

注意事项 (Note)

For research use only .