

产品名称: NUP62CL Rabbit Monoclonal Antibody
产品货号: AMRe02826



产品概述 (Summary)

产品名称 (Production Name)	NUP62CL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NUP62CL
别名 (Alternative Names)	NUP62L; Nucleoporin-62 C-terminal-like protein
基因 ID (Gene ID)	54830
蛋白 ID (SwissProt ID)	Q9H1M0.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 21 kDa

研究背景 (Background)

This gene encodes a protein containing a domain found in nucleoporins which are glycoproteins found in nuclear pore


产品名称: NUP62CL Rabbit Monoclonal Antibody
产品货号: AMRe02826



complexes. Alternative splicing results in multiple transcript variants.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NUP62CL in HeLa, Ramos lysates using NUP62CL antibody.

注意事项 (Note)

For research use only .