

产品名称: E1B AP5 Rabbit Monoclonal Antibody
产品货号: AMRe02828



产品概述 (Summary)

产品名称 (Production Name)	E1B AP5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HNRNPUL1
别名 (Alternative Names)	E1BAP5; E1B-AP5; HNRPUL1
基因 ID (Gene ID)	11100
蛋白 ID (SwissProt ID)	Q9BUJ2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 96 kDa; Observed MW: 118 kDa

研究背景 (Background)

Acts as a basic transcriptional regulator. Represses basic transcription driven by several virus and cellular promoters. When

产品名称: E1B AP5 Rabbit Monoclonal Antibody
产品货号: AMRe02828




associated with BRD7, activates transcription of glucocorticoid-responsive promoter in the absence of ligand-stimulation. Plays also a role in mRNA processing and transport. Binds avidly to poly(G) and poly(C) RNA homopolymers in vitro.


研究领域 (Research Area)

Microbiology


图片 (Image Data)



Immunocytochemistry analysis of E1B AP5 (green) in Jurkat using E1B AP5 antibody, and DAPI (blue).



Western blot analysis of E1BAP5 in A549, HL-60, U251 lysates using E1BAP5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using E1BAP5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: E1B AP5 Rabbit Monoclonal Antibody
产品货号: AMRe02828



注意事项 (Note)

For research use only .